CIRUGIA GUIADA A SU MANERA

exocad

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4.6 mm Master Cylinder H6

Article Number	CGS-QMC-10
Inner Diameter ϕ	5.11 mm
Height H	5.99 mm
To platform: 10.53 mm	
To bone level: 10.53 mm	

Biocorizons® - Guided Surgery Kit CGS4000

1	2	3	4	5	6	7
ϕ 4.6 mm	ϕ 2 mm L 24 mm	ϕ 2.5 mm L 24 mm	ϕ 3.2 mm L 24 mm	ϕ 3.7 mm L 24 mm	ϕ 4.1 mm L 24 mm	Stop Position SP3

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Tools and Guides:

- Mucosa Punch, ϕ 4.0 mm, guided, (034.011)
- Drill Handle, self-locking, ϕ = 3.5 mm, h = 3 mm (034.294)
- Drill Handle, self-locking, ϕ = 2.2 mm, h = 3 mm (034.291)
- Drill Handle, self-locking, ϕ = 3.5 mm, h = 3 mm (034.294)
- Drill Handle, self-locking, ϕ = 3.7 mm, h = 3 mm (034.295)
- Milling Cutter, ϕ 3.5 mm, guided (034.415)
- X Pilot VeloDrill™, guided, ϕ 2.2 mm, L 20 mm (066.1501)
- X VeloDrill™, guided, ϕ 3.5 mm, L 20 mm (066.1504)
- X VeloDrill™, guided, ϕ 3.7 mm, L 16 mm (066.1305)
- BLX Guided Implant Driver, ratchet, SST, (066.4404)

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Implant systems (83)

35 **straumanngroup**

- BLX - Bone Level X Roxolid® SLActive®
- ϕ 4.5 mm RB, 10 mm
- Guided BLX Sleeves Fully Guide... For Hard Bone Density with: Str...
- Select prosthetic component

32 **straumanngroup**

- BLX - Bone Level X Roxolid® SLActive®
- ϕ 3.75 mm RB, 12 mm
- Select surgical component
- Select prosthetic component

42 **Tooth 42**

- Select implant

Manufacturer: Straumann® Group

System: BLX - Bone Level X Roxolid® SLActive®

Implant Analog: BLX, ϕ 3.75 mm RB, SLActive® 12 mm, Roxolid®

Article No.: 061.4312

Length: 12 mm

Platform ϕ : 3.5 mm

Body ϕ : 3.75 mm

Comment: Material: Roxolid®, Surface: SLActive®

Buttons: SLEEVES / KITS > PROSTHETICS > CONFIRM

Message: You can now place the implant with the mouse.

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DATA VIEW SURFACE GENERATION

Study: Dual-Scan-Protocol, Example

Series: Patient with markers in both jaws (407 Images)

Modality: CT
Description: Patient with markers in both jaws
No. of images: 407
Image type: AXIAL
Series date: 5/15/2021

Select file set Load series Discard series

Incognito mode

Approval of Planning

The implant planning is now completed. You are now kindly asked to approve your planning and confirm the items below before the implant planning report and data files are generated.

Patient Information		
Name	Project	DICOM
Birth date		


- I agree to study the implant planning report before commencing surgery and verify the provided information with regard to the surgery.
- I confirm that I have performed the planning with due care and that I am satisfied with the medical and clinical results of the planning. I chose the most appropriate implant solution and observed safety margins to any existing dentition, existing and intended dental restorations, and anatomical structures.

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Screenshot and Image Management

Screenshots & Images

+ Add all images to report Save all images Delete all images

Caption: Implant cross 1 (Anchor Pin 1) Included in report: (1) <input checked="" type="checkbox"/> Implant planning report <input type="checkbox"/> Surgical report	Caption: Panorama Included in report: (1) <input type="checkbox"/> Implant planning report <input checked="" type="checkbox"/> Surgical report	Caption: Implant cross 1 (Tooth 44) Included in report: (1) <input checked="" type="checkbox"/> Implant planning report <input type="checkbox"/> Surgical report	Caption: Implant cross 1 (Tooth 34) Included in report: (1) <input checked="" type="checkbox"/> Implant planning report <input type="checkbox"/> Surgical report

Add screenshot or load image

Capture main view only Capture manually Capture all views Load from file

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