

exocad GmbH

Natalia Gonsior
Senior Marketing Specialist
Rosa-Parks-Straße 2
64295 Darmstadt, Germany
Phone +49 (0)6151 6294132
ngonsior@exocad.com
[exocad.com](https://www.exocad.com)

XSNAP IS EXOCAD'S NEW ADD-ON MODULE FOR PRINTABLE ARTICULATORS

Perform natural jaw movements with a 3D-printed jaw model

DARMSTADT, Germany – February 8, 2023 – exocad, an Align Technology, Inc. company and a leading dental CAD/CAM software provider, now provides dental models with a 3D-printed articulator. The new xSNAP add-on module for *DentalCAD 3.1 Rijeka* allows users to simulate physiological mandibular movements (protrusion, mediotrusion and laterotrusion) without the use of a conventional articulator. xSNAP is also the first 3D joint system for dental models produced in a 3D printer, which enables the simulation of dynamic jaw movements. The xSNAP attachments are equipped with integrated articular paths and condylar path inclinations of 15, 25, 35 and 45 degrees. Additionally, users can choose between a flat incisal plate and an incisal plate with a 30-degree inclination. exocad offers added value to many users with the easy click mechanism for joining the models and precisely designed guideways.

“We’re happy to offer an even more cost-effective and accelerated workflow for dental models to further empower our customers,” said Novica Savic, CCO and general manager at exocad.

exocad users also benefit from these additional software features:

- The co-printed component eliminates the need to subsequently articulate acrylic models
- Joint guidance is dynamically adjustable for each printer and personal preferences
- No additional support structures are necessary because of the adaptive support of the xSNAP attachment

The *xSNAP Module* is available as a library for exocad’s *Model Creator* module in *DentalCAD 3.1 Rijeka*, the newest generation of exocad’s open architecture CAD software for labs and

Press Release

clinics. Additional information is available at exocad.com/our-products/exocad-dentalcad/xsnap-module

###

About exocad

exocad GmbH, an Align Technology, Inc. company, is a leading dental CAD/CAM software provider. exocad vigorously pushes the boundaries of digital dentistry, providing flexible, reliable, and easy-to-use CAD/CAM software for dental labs and dental practices worldwide. More than 50,000 valued customers plan implants and create functional and refined restorations with exocad's *DentalCAD*, *ChairsideCAD* and *exoplan* software. exocad and *DentalCAD*, among others, are trademarks of exocad GmbH or one of its subsidiaries or affiliated companies, and are registered in the U.S. and other countries. For more information and a list of exocad reseller partners, please visit exocad.com.